TOP NEWS

05

2016

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***FULLY AUTOMATED FORMALIN MIXING STATION –
THE SOLUTION FOR CONTACT-FREE FORMALIN
MIXING***

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Various institutions have published information on safe handling of formaldehyde and formalin since the reclassification of formaldehyde to carcinogenic and germ cell mutagenic by the European Union. The institutions advice, among other things, to avoid the mixture of a 3.7% formalin solution of 37% stock solution and instead to convert to ready for use 3.7% formalin solution. This advice however, cannot be implemented in practice by all pathologies. We have developed a fully automated formalin mixing station to tackle that issue. Hence, lab staff does not get into contact with 37% stock solution and the formalin concentration in the pathology is subsequently reduced.

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The **automatic formalin mixing station** is designed for **contact-free, easy and reliable preparation of formalin solutions** in laboratories. Thanks to its **compact design**, the mixing station is also suitable for small laboratories. The device is to be connected directly to the water supply. It prepares the mixture within the tank in the desired mixing ratio. The integrated microprocessor within the device allows you to easily control all operations, such as **preparing a 3.7% standard solution**, through one single surface at the control panel with LCD display. Additionally, it is possible to prepare a **customized mixture by using variable quantities of formalin, purified water and a buffer liquid**.

FULLY AUTOMATED FORMALIN MIXING STATION

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During the **automatic mixing process**, 37% stock solution is being pumped from the front rack into the tank, then water and a buffer solution. These liquids are subsequently **mixed for 15 minutes**. Additionally, the **remaining solution is blended automatically 6 times a day for 15 minutes. Electronic and mechanical safety precautions** keep liquids from overflowing within the tank.

The prepared mixture can be drained directly from a tap at the device or transferred to multiple remote stations that are up to 200 meters away.

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EUROPEAN MAC VALUES FOR FORMALDEHYDE

As a result of the reclassification of Formaldehyde to carcinogenic (cancer-causing) and germ cell mutagenic by the European Union, the formally known MAC value for working involving Formaldehyde, $0,3 \text{ ml/m}^3$ or $0,37 \text{ mg/m}^3$ and the short term value for short term peak contamination of $0,74 \text{ mg/m}^3$ shall not be exceeded by any means in Germany.

Please inform yourself as soon as possible about the maximum concentration values in your country.

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Operation: control panel with LCD display
Mains supply: 230 Volt 10 Amp / 50-60 Hz
Mains connection: shockproof plug
Power consumption: 250 W
Water connection: DN 15
Usable liquids: formalin, buffer liquid, distiller water or drinking water
Water overflow detecting sensor
Audible notification
Notification LEDs: red and green
Minimum liquid volume for standard solution: 10 liters
Minimum liquid volume for customized solution: 1 liter/liquid
Liquid level monitoring precision: 0.1 liter (+/- 1%)
Minimum residue liquid level required in the tank: 5 liters
Tank volume: 100 liters
Security tank volume: 67 liters
Dimensions W/L/H: 60 x 70 x 1320 mm
Dead weight: 100 kg
Ambient temperature: 5°C to 45°C
Type: AFMDS-100
Article No.: 099.900.196

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Please refer to our brochure and leaflet for more information about the reclassification of formaldehyde and our exhaust air units with which you are able to comply with the maximum concentration values at the workplace for work involving formaldehyde and formalin!

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Straße 2 A
DE – 93049 Regensburg
GERMANY

Phone +49 (0) 9 41/20 86 48-0
Fax +49 (0) 9 41/20 86 48-29
E-Mail info@KUGEL-medical.de
www.KUGEL-medical.de

KUGEL
medical
■ ■ ■ ■