03/2018

TOP NEWS

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***CLEANCUBE CLASSIC: YOUR ALLROUND TALENT FOR
SURFACE DISINFECTION***

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Break the chain of infection through modern 3D disinfection without risk to health!

Imagine a dry fog on all surfaces with no harmful residues. The Swiss Company CleanService Group AG developed a new system for large surface disinfection with its CleanCube. The hunt for hidden spots, the struggle to cover all surfaces, is history using 3D disinfection. Your employees are free to cover tasks that are more challenging.

A fog of ultrafine droplets covers all accessible surfaces within a room thanks to the constant stream of airflow. Nothing is left behind, not even where access is most difficult. As the particles of the fog are extremely tiny, they increase substantially the surface contact, attacking bacteria membranes exponentially. This 3D disinfection process will normally take just a few minutes for a standard sized room. After a short waiting time to allow all particles in the air to settle and to impact on the micro-organisms inside the room, the latter can be safely entered again, no rinsing or post-cleaning is necessary.

CSGs 3D disinfection proves to be successful in fighting against resistant micro-organisms such as MRSA and nosocomial infections, thus achieving a reliable, flawless, time- and cost-saving air and surface decontamination.

CLEANCUBE CLASSIC: YOUR ALLROUND TALENT FOR SURFACE DISINFECTION

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The versatile all-rounder at a glance:

- Precise Reproducibility with pump driven dosing
- Only 32 x 32 x 36 cm
- Large room size up to 700 m³
- Status visible from far with the multicolor ring light
- Programmable stepping and timer function
- Optional extension tube, flange and gun

Working Principle:

The active substance used is **hydrogen peroxide**, an **environment friendly substance**. In a complex manufacturing procedure, the active substance is stabilized and **boosted with silver**, thus achieving a greatly improved effectiveness against micro-organisms. The traces of silver remaining on the treated surfaces are **not visible and non toxic**. However, they efficiently inhibit a renewed contamination.

The elemental oxygen (O₂) separated by the hydrogen peroxide **attacks the cell walls of the micro-organisms directly**. The chemical reaction between the oxygen and the cell wall molecules will cause these to be denatured and destroyed.

This effect is intensified by the silver ions which form a bond with the disulfide bridge of certain proteins of micro-organisms, thereby **inactivating or precipitating these proteins**.

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CleanLine 015 has been proven effective against:

- gram positive / gram negative bacteria
- spore-forming organisms
- viruses (with and without protective coat)
- Yeast
- a wide spectrum of protozoa

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