

Preparation Cabinet NSP-1000 – NSP-1500 G with Air Suction System



The sturdy construction is completely made of brushed stainless steel (material no. 1.4301 - ANSI 304). The cabinet has a 100 mm high surrounding base cover and 2 double-walled stainless steel doors with 2 big bow handles. 5 shelf inserts can be placed inside the cabinets and height-adjusted by means of grid rails. Please note that the exhaust air system must always be plugged in so that enough air can flow through to discharge safely escaping formalin vapors.

Technical Data:

Dimensions W/D/H:	1000 x 600 x 2100 mm
Exhaust Air capacity:	150 m ³ /h
Material:	Stainless steel (Material No. 1.4301 - ANSI 304)
Doors:	2 doors made of stainless steel
Exhaust duct diameter:	150 mm
Type:	NSP-1000
Article No.:	090.020.100

Additional models see back page

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg / Germany

Phone +49/941/208648-0
Fax +49/941/208648-29
www.KUGEL-medical.de



KUGEL
medical
■ ■ ■ ■

Preparation cabinet NSP-1000 – NSP-1500 G with Air Suction System

Model NSP-1200:

Dimensions W/D/H:	1200 x 600 x 2100 mm
Exhaust Air capacity:	Min. 200 m ³ /h
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Doors:	2 doors made of stainless steel
Exhaust duct diameter:	150 mm
Type:	NSP-1200
Article No.:	090.020.110

Model NSP-1500:

Dimensions W/D/H:	1500 x 600 x 2100 mm
Exhaust Air capacity:	Min. 250 m ³ /h
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Doors:	2 doors made of stainless steel
Exhaust duct diameter:	150 mm
Type:	NSP-1500
Article No.:	090.020.120

Model NSP-1000-G:

Dimensions W/D/H:	1000 x 600 x 2100 mm
Exhaust Air capacity:	Min. 150 m ³ /h
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Doors:	2 doors made of glass
Exhaust duct diameter:	150 mm
Type:	NSP-1000-G
Article No.:	090.020.105

Model NSP-1200-G:

Dimensions W/D/H:	1200 x 600 x 2100 mm
Exhaust Air capacity:	Min. 200 m ³ /h
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Doors:	2 doors made of glass
Exhaust duct diameter:	150 mm
Type:	NSP-1200-G
Article No.:	090.020.115

Model NSP-1500-G:

Dimensions W/D/H:	1500 x 600 x 2100 mm
Exhaust Air capacity:	Min. 250 m ³ /h
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Doors:	2 doors made of glass
Exhaust duct diameter:	150 mm
Type:	NSP-1500-G
Article No.:	090.020.125