

Autopsy Table ST 10/200 – Rotating



Model ST 10/100 with additional tap beside the basin

The autopsy table is mounted on a center support column and can be rotated to the right and left by approx. 180°. The autopsy table is completely made of stainless steel (material no. 1.4301- ANSI 304). The work surface features a high, profiled surround and a seamlessly welded basin, 400 x 500 x 200 mm, is with a tube overflow and a sieve insert is installed at one end of the table. Additionally, the table top is equipped with a waterproof sound-proofing and features a negative inclination arranged diagonally to the drain. All power connections are built into the base and can be easily accessed by means of an inspection opening.

The Autopsy Table ST 10/150 model includes the following features:

- Knee-operated mixing tap for cold and warm water at the center column
- 3 m shower hose with hand sprinkler
- 2 splash proof electrical outlets
- Electrical and water connections are preinstalled
- UP / DOWN button for height adjustment by means of an integrated electro mechanical drive
- The safety requirements for drinking water as per DIN 1988 part 4 and EN DIN 1717 are to be provided by the customer

Technical Data:

Dimensions W/D/H:	2600 x 850 x 850 mm
Voltage:	230 Volt / 1 Ph / 50 Hz
Carrying capacity:	max. 300 kg
Type:	ST 10/200
Article No.:	010.010.200

Alternatively available with a carrying capacity of 450 kg at an extra charge – Article No.: 010.010.858

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

