

# Autopsy Table ST 10/350

## Fixed Mounted – BASIC line



The table is set on an installation column and two height-adjustable feet to equalize out any unevenness. The autopsy table is completely made of stainless steel (material no. 1.4301- ANSI 304). The work surface features a high, profiled surround and a seamlessly welded basin, 400 x 500 x 200, with a tube overflow are installed to the end of the table. The table top has a negative inclination arranged diagonally to the drain. The drains from the supporting surface and the basin are brought together in one drain connection. All power connections are built into the base and can be easily accessed by means of an inspection opening.

### **The Autopsy table ST 10/350 model includes the following features:**

- Knee-operated mixing tap for cold and warm water
- 3 m shower hose with hand sprinkler
- 2 splash proof electrical outlets 230 Volt - German style
- The safety requirements for drinking water as per DIN 1988 part 4 and EN DIN 1717 are to be provided by the customer
- Electrical and water connections inside the installation column are to be provided by the customer

### **Technical Data:**

Dimensions W/D/H:	2600 x 850 x 850 mm
Voltage:	230 Volt / 1 Ph / 50 Hz
Carring capacity:	max. 300 kg
Type:	ST 10/350
Article No.:	010.010.350

**KUGEL medical GmbH & Co. KG**  
Hermann-Köhl-Str. 2A  
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0  
Fax +49/941/20 86 48-29  
[www.KUGEL-medical.de](http://www.KUGEL-medical.de)



**KUGEL**  
medical  
■ ■ ■ ■

# Autopsy Table ST 10/350

## Fixed Mounted – BASIC line



Removable stainless steel box  
Dimensions: W/D/H 830 x 200 x 60 mm  
Type: AB-ST 350-830  
Article No. 010.010.355



Stainless steel body supports (4 pc.)  
Dimensions: W/D/H 830 x 200 x 60 mm  
Type: KAB ST-10-30  
Article No. 010.021.060

Neck rest made of black rubber  
Dimensions: W/D/H 150 x 120 x 120 mm  
Type: NS 10  
Article No. 010.020.080

**KUGEL medical GmbH & Co. KG**  
Hermann-Köhl-Str. 2A  
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0  
Fax +49/941/20 86 48-29  
[www.KUGEL-medical.de](http://www.KUGEL-medical.de)



**KUGEL**  
medical  
■ ■ ■ ■