

Autopsy Table ST 10/500 – Moveable



The autopsy table is made completely of stainless steel (material no. 1.4301). The substructure consists of rugged, torsion-resistant stainless steel profiles with 4 steerable solid rubber wheels (D= 125 mm), with permanently lubricated, maintenance-free special roller bearings, 2 of them lockable. The supporting surface features a high, profiled surround and a seamlessly integrated basin, 400 x 500 x 200, with overflow pipe and sieve insert. The table top is equipped with sound-proofing and features a negative inclination arranged diagonally to the drain. The drains from the supporting surface and the basin are brought together in one drain connection. Delivery includes a siphon.

Technical Data:

Dimensions W/D/H:	2600 x 750 x 850 mm
Carrying capacity:	max. 200 kg
Type:	ST 10/500
Article No.:	010.010.500

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

