

Bench vice - S



The entire bench vice is made of stainless steel (material no. 1.4301 – ANSI 304) with prismatic precision channel and a sturdy shaft with acme screw thread for quick open and close. Round clamping jaws with 6 teeth (not replaceable). Delivery including wall mount made of stainless steel.

Technical details:

Diameter of jaw:	approx. 120 mm
Clamping area:	approx.. 0 – 390 mm
Dimensions W/D/H:	607 x 542 x 930 mm
Type:	SST-200x120-400
Article No.	019.300.027

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



KUGEL
medical
■ ■ ■ ■