

Cabinet for acid and alkalis with pull-out shelves and PP trays



Image: cabinet for acid and alkalis with 3 pull-out trays

The cabinet for acids and alkalis inclusive 4 or 6 pull-out shelves and PP trays is particularly suitable for a safe storage of corrosive hazardous materials, especially acids and alkalis according to DIN EN 14727. It is not suitable for the storage of flammable liquids. The Outer corpus is made of sheet steel.

The pull-out-shelves with trays have a sturdy construction. Each shelf is equipped with a removable PP tray, impervious to liquids and highly resistant against corrosives. Maximum shelf load: 25 kg (evenly distributed).

Ventilation of both compartments through a single spigot (diameter 75 mm). The cabinet can be integrated into existing master-key systems.

Technical Data:

Outside Dimensions W/D/H:	approx. 600 x 605 x 1970 mm
Dimensions of each Tray B/D/H:	approx. 485 x 565 x 930 mm (4 compartments)
Dead weight:	approx. 125 kg
Carrying weight drawer:	25 kg
According to:	DIN EN 14727
Exhaust duct diameter:	75 mm
Colour:	Lightgrey RAL 7035
Type:	SL.196.60-4
Article No.:	090.020.585

Additional models see back page

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



Cabinet for acid and alkalis with pull-out shelves and PP trays

Model SL.196.60-6:

Outside Dimensions W/D/H:	approx. 600 x 605 x 1970 mm
Dimensions of each Tray B/D/H:	approx. 485 x 565 x 930 mm (6 compartments)
Dead weight:	approx. 125 kg
Carrying weight drawer:	25 kg
According to:	DIN EN 14727
Exhaust duct diameter:	75 mm
Colour:	Lightgrey RAL 7035
Type:	SL.196.60-6
Article No.:	090.020.590

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

