


Corian cutting board



The cutting board is completely made of corian. The top thickness is approx. 12 mm. The board includes 4 small feet approx. 25 mm height and 20 mm diameter. The cutting board is not suitable for thermal sterilization.

Technical Details:

Dimensions W/D:	500 x 400 mm
Material:	12 mm Corian
Application:	black point with 140 mm diameter
Groove:	ca. 16 mm wide and 10 mm low
Type:	SE-CO-50/40-SE-R
Article-No.:	060.012.002

Accessoires:

Ruler made of stainless steel with cm and mm measurements	
Length:	300 mm
Material:	Stainless Steel
Type:	LIN-SB 300
Article-No.:	060.019.031

Additional models see back page

Special bodies are available on inquiry

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2a
DE - 93049 Regensburg

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



Corian cutting board

Additional models

Corian cutting board

Dimensions W/D: 500 x 400 mm
Material: 12 mm Corian
Typ: SE-CO-50/40
Article-No.: 060.012.009

Corian cutting board with groove

Dimension W/D: 500 x 400 mm
Material: 12 mm Corian
Groove: ca. 16 mm wide and 10 mm low
Typ: SE-CO-50/40-R
Article-No.: 060.012.003

Corian cutting board with black point

Dimension W/D: 500 x 400 mm
Material: 12 mm Corian
Application: round point with 140 mm diameter
Typ: SE-CO-50/40-SE
Article-No.: 060.012.012

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2a
DE - 93049 Regensburg

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

