

Polyethylene cutting board - blue



The cutting board is completely made of polyethylene. The top thickness is approx. 20 mm. The board includes 4 small feet approx. 25 mm height and 15 mm diameter. not suitable for thermal sterilization

Technical Details:

Dimensions W/D:	500 x 400 mm
Material:	20 mm PE-Material
Colour:	blue
Type:	SB-PE-50/40-BL
Article-No.:	060.012.119

Accessoires:

Ruler made of stainless steel with cm and mm measurements

Length:	300 mm
Material:	Stainless Steel
Type:	LIN-SB 300
Article-No.:	060.019.031

Special dimension are available – please ask for more information

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2a
DE - 93049 Regensburg

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



KUGEL
medical
■ ■ ■ ■