


Polyethylene cutting board



The cutting board is completely made of polyethylene. The top thickness is approx. 20 mm. The board includes 4 small feet approx. 15 mm height and 25 mm diameter. The cutting board is not suitable for thermal sterilization.

Technical Details:

Dimensions W/D:	500 x 400 mm
Material:	20 mm PE-Material
Type:	SE-PE-50/40
Article-No.:	060.012.010

Accessoires:

Ruler made of stainless steel with cm and mm measurements

Lenght:	300 mm
Material:	Stainless Steel
Type:	LIN-SB 300
Article-No.:	060.019.031

Additional models see back page

Special bodies are available on inquiry

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2a
DE - 93049 Regensburg

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



Polyethylene cutting board

Polyethylene cutting board with groove

Dimensions W/D:	500 x 400 mm
Material:	20 mm PE-Material
Groove:	ca. 16 mm wide and 10 mm low
Type:	SE-PE-50/40-R
Article-No.:	060.012.013

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2a
DE - 93049 Regensburg

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

