

# Digital Scales with Flexible Display



Scale type 7755



Scale type 7755 Plus with fixed mounted holding device

Completely electronic digital scale for quick, reliable and easy weighing. Taring over the entire weighing range with only one button. Sturdy housing made of grey painted steel and a weighing platform made of stainless steel with flexible display made of sturdy plastic. The scale can be either operated with the integrated battery or by connection if to a 230 Volt 50/60 Hz power supply. Working hours of the battery approximately 8 hours. The weighing platform and the display are connected with a coiled cord, length approximately 1,5 m.

## **Standard delivery includes:**

- Type 7755 / 7756 / 7757 with flexible display incl. small holding device and power pack
- Type 7755-PLUS with fixed mounted holding device for the display and power pack

## **Technical Data:**

|                                   |                   |
|-----------------------------------|-------------------|
| Dimensions W/D/H:                 | 310 x 275 x 35 mm |
| Dimensions display housing W/D/H: | 260 x 140 x 70 mm |
| Weighting range:                  | 0 - 15 kg         |
| Accuracy:                         | 5 gram            |
| Weight:                           | 4 kg              |
| Type                              | 7755-15           |
| Article No.:                      | 090.010.065       |

*Additional models see back page*

**KUGEL medical GmbH & Co. KG**  
Hermann-Köhl-Str. 2A  
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0  
Fax +49/941/20 86 48-29  
[www.KUGEL-medical.de](http://www.KUGEL-medical.de)



# Digital Scales with Flexible Display

## **Model 7755-30:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 310 x 275 x 35 mm |
| Weighting range:  | 0 - 30 kg         |
| Accuracy:         | 10 gram           |
| Weight:           | 4 kg              |
| Type:             | 7755-30           |
| Article No.:      | 090.010.070       |

## **Model 7755-60:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 310 x 275 x 35 mm |
| Weighting range:  | 0 - 60 kg         |
| Accuracy:         | 20 gram           |
| Weight:           | 4 kg              |
| Type:             | 7755-60           |
| Article No.:      | 090.010.065       |

## **Modell7755-PLUS-15:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 310 x 275 x 35 mm |
| Weighting range:  | 0 - 15 kg         |
| Accuracy:         | 5 gram            |
| Weight:           | 6 kg              |
| Type:             | 7755-PLUS-15      |
| Article No.:      | 090.010.079       |

## **Model 7755-PLUS-30:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 310 x 275 x 35 mm |
| Weighting range:  | 0 - 30 kg         |
| Accuracy:         | 10 gram           |
| Weight:           | 6 kg              |
| Type:             | 7755-PLUS-30      |
| Article No.:      | 090.010.080       |

## **Model 7755-PLUS-60:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 310 x 275 x 35 mm |
| Weighting range:  | 0 - 60 kg         |
| Accuracy:         | 20 gram           |
| Weight:           | 6 kg              |
| Type:             | 7755-PLUS-60      |
| Article No.:      | 090.010.081       |

**KUGEL medical GmbH & Co. KG**  
Hermann-Köhl-Str. 2A  
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0  
Fax +49/941/20 86 48-29  
[www.KUGEL-medical.de](http://www.KUGEL-medical.de)



# Digital Scales with Flexible Display

## **Model 7756-60:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 520 x 400 x 70 mm |
| Weighting range:  | 0 - 60 kg         |
| Accuracy:         | 20 gram           |
| Weight:           | 16 kg             |
| Type:             | 7756-60           |
| Article No.:      | 090.010.085       |

## **Model 7756-150:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 520 x 400 x 70 mm |
| Weighting range:  | 0 - 150 kg        |
| Accuracy:         | 50 gram           |
| Weight:           | 16 kg             |
| Type:             | 7756-150          |
| Article No.:      | 090.010.086       |

## **Model 7757-150:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 660 x 545 x 85 mm |
| Weighting range:  | 0 - 150 kg        |
| Accuracy:         | 100 gram          |
| Weight:           | 28 kg             |
| Type:             | 7757-150          |
| Article No.:      | 090.010.090       |

## **Model 7755-PLUS-30:**

|                   |                   |
|-------------------|-------------------|
| Dimensions W/D/H: | 660 x 545 x 85 mm |
| Weighting range:  | 0 - 200 kg        |
| Accuracy:         | 200 gram          |
| Weight:           | 28 kg             |
| Type:             | 7757-200          |
| Article No.:      | 090.010.091       |

**KUGEL medical GmbH & Co. KG**  
Hermann-Köhl-Str. 2A  
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0  
Fax +49/941/20 86 48-29  
[www.KUGEL-medical.de](http://www.KUGEL-medical.de)

