

Hand eye shower HAD 020



Broad water stream shower head inclined by 45°, low-scaling spray plate made of plastic
Rubber protection with dust lid and a eye shower symbol
Does not turn off naturally – comprises a lockable deployment handle with integrated flow regulator
8L/min

Technical Data:

Hose with stainless steel sheath:	1,5 m
Table duct, dimensions:	28 x 1,5 x 80 mm
Fluid-carrying parts:	made of brass
Handle and deployment handle:	made of plastic
Incl. eye shower sign, self-adhesive 100 x 100 mm	
Water flow volume:	approximately 8L/min
Water supply:	DN 15
Water pressure:	min. 2 bar
Type:	HAD 020
Article No.:	060.200.020

Standards and Guidelines:

The eye shower meets the EN 15154, part II
Safety markings according to ISO 3864-1

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE – 93049 Regensburg / Germany

Phone +49/941/208648-0
Fax +49/941/208648-29
www.KUGEL-medical.de

