

Electronic Hanging Scale Construction

EHW-20-ASG



A very sturdy hanging scale for minor loads with removable weighing bowl made of stainless steel. The scale is being delivered with batteries (3 x 1,5 V). for immediate use. An AUTO-OFF function is being installed which automatically turn off the scale if there is no change of load after 3 minutes to preserve the battery. Delivery includes frame construction completely made of stainless steel with 4 rubber feet. This scale is perfect in conjunction with an autopsy table.

Technical Data:

Dimensions W/D/H:	550 x 425 x 850mm
Weighing range:	max. 15 kg
Accuracy:	20 gram
Type:	EHW-20-ASG
Article-No.:	090.012.045

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

