

Instrument Trolley BT 80/6060 - 9560



The instrument trolley is completely made of stainless steel (Material no. 1.4301 – ANSI 304). The supporting surface features a high, profiled surround. The trolley is equipped with 4 steerable wheels.

Technical Data:

Dimensions W/D/H: 600 x 600 x 800 mm
Wheels: Ø 100 mm
Caring capacity: max. 100 kg
Material: Stainless steel (Material No. 1.4301 – ANSI 304)
Type: BT 80/6060
Article no.: 080.050.010

Additional models see back page

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de



Instrument Trolley BT 80/6060 - 9560

Modell BT 80/7560:

Dimensions W/D/H:	750 x 600 x 800 mm
Wheels:	Ø 100 mm
Caring capacity:	max. 100 kg
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Type:	BT 80/75060
Article no.:	080.050.020

Modell BT 80/9560:

Dimensions W/D/H:	950 x 600 x 800 mm
Wheels:	Ø 100 mm
Caring capacity:	max. 100 kg
Material:	Stainless steel (Material No. 1.4301 – ANSI 304)
Type:	BT 80/9560
Article no.:	080.050.030

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

