




Portable waterproof autopsy saw SF-4000-P

Tragbare wasserdichte Autopsiesäge SF-4000-P



Pros

- ⊕ efficient 13200 oscillation per minute cycle
- ⊕ unique ergonomic easy grip design
- ⊕ watertight to IP67 rating, easy to clean and decontaminate
- ⊕ operates on safe 40 Volts
- ⊕ light weight 1300g hand piece
- ⊕ comprehensive chrome plated blade range
- ⊕ available in 110V and 240V versions
- ⊕ with 4.0 m power cord
- ⊕ anodized aluminium body with 5 bearing sealed mechanism
- ⊕ vibration value 4.12 m/s²
- ⊕ 37 db(A) noise level (no blade)
- ⊕ typical load 1.5 Amp
- ⊕ DC, permanent magnet motor
- ⊕ transformer manufactured to IEC – EN60742
- ⊕ enclosure manufactured to IEC529 and DIN40050

Vorteile

- ⊕ effiziente Oszillation: 13200 Zyklen pro Minute
- ⊕ einzigartiges ergonomisches Griffdesign
- ⊕ wasserfest gemäß IP67, einfach zu reinigen und desinfizieren
- ⊕ sicherer 40 Volt Betrieb
- ⊕ leichtes 1300g Handstück
- ⊕ breites Sortiment an verchromten Sägeblättern
- ⊕ verfügbar als 110V und 240V Version
- ⊕ mit 4,0 m Stromkabel
- ⊕ eloxiertes Aluminiumgehäuse mit 5 Lager abgedichtetem Mechanismus
- ⊕ Vibrationswert 6,93 m/s²
- ⊕ Lärmpegel 37 db (A) (ohne Sägeblatt)
- ⊕ typische Belastung 1,5A
- ⊕ Motor mit Gleichstrom Permanentmagnet
- ⊕ Umspanner hergestellt nach IEC – EN60742
- ⊕ Gehäuse hergestellt nach IEC529 und DIN40050