

Safety Tall Cabinet with one Wing Door



The fire-resistant wing door cabinet according to DIN EN 14470-1 (Type 90), DIN 12925-1 (FWF90) and TRbF 20 (attachment L), especially designed for the safe and correct storage of dangerous liquids in work areas. Complete closing mechanics as well as all safety relevant construction units are installed outside of the cabinet interior. Grounding connection of each drawer. The cabinet includes 3 drawers.

Technical Data:

Outside Dimensions W/D/H:	600 x 615 x 1298 mm
Inside Dimensions W/D/H:	450 x 520 x 1070 mm
Absorb volume:	approx. ca. 5 ltr.
Dead weight:	approx. 170 kg
Carrying weight drawer:	each 25 kg
According to:	DIN EN 14470-1
Exhaust duct diameter:	75 mm
Air flow:	approx. 2 m ³ / hour
Colour:	Lightgrey RAL 7035
Type:	SIHS-600
Article No.:	090.020.520

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

