

Safety Tall Cabinet with two Wing Doors



The fire-resistant wing door cabinet according to TIN EN 14470-1 (Type 90), DIN 12925-1 (FWF90) and TRbF 20 (attachment L), especially designed for the safe and correct storage of dangerous goods in work areas. Complete closing mechanics as well as all safety relevant construction units are installed outside of the cabinet interior.

Technical Data:

Outside Dimensions W/D/H:	1200 x 615 x 1968 mm
Inside Dimensions W/D/H:	1050 x 520 x 1740 mm
Dead weight:	approx. 420 kg
According to:	DIN EN 14470-1
Exhaust duct diameter:	75 mm
Air flow:	approx. 9 m ³ / hour
Colour:	Lightgrey RAL 7035
Type:	SIHS-1200
Article No.:	090.020.500

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

