

Safety Cabinet with one Drawer



The fire-resistant door cabinet according to DIN EN 14470-1 (Type 90), DIN 12925-1 (FWF90) and TRbF 20 (attachment L), especially designed for the safe and correct storage of dangerous liquids in work areas. Fully extendable, low-friction and lockable in each position. The drawer is closing by thermal release in case of a fire and has a grounding connection on the outside.

Technical Data:

Outside Dimensions W/D/H:	1100 x 570 x 630 mm
Inside Dimensions W/D/H:	980 x 450 x 500 mm
Useable inside height:	approx. 460 mm
Absorb volume:	approx. 27 ltr.
Dead weight:	approx. 165 kg
Carrying weight drawer:	50 kg
According to:	DIN EN 14470-1
Exhaust duct diameter:	75 mm
Air flow:	approx. 2m ³ / hour
Colour:	Lightgrey RAL 7035
Type:	SIHS-1100
Article No.:	090.020.420

KUGEL medical GmbH & Co. KG
Hermann-Köhl-Str. 2A
DE - 93049 Regensburg/Germany

Phone +49/941/20 86 48-0
Fax +49/941/20 86 48-29
www.KUGEL-medical.de

