





# Grossing Station with Integrated Laminar Airflow Supply Element **WMT 2250-E-BTR Model**

FP001068 | 060.300.210



Wall-mounted grossing station equipped with a laminar airflow ventilation system to provide optimal operator protection against exposure to chemical vapours during the handling and examination of tissue specimens.



## Technical specifications

- Dimensions: **2250 x 850 x 2100 mm**
- Electrical height adjustment : **780- 1060 mm**
- Tabletop in **stainless steel 316** | Structure in **stainless steel 304**
- **Integrated laminar airflow supply element (ductable)**
- With integrated dimmable **4000K LED lighting**
- **Peripheral ventilation** around the workspace (exhaust duct located at the center of the extraction area)
- **Rear extraction** behind the basin
- Work surface consisting of **one removable, non-perforated plate**
- Seamless welded basin **400 x 400 x 200 mm x 3 ½ "** on the right side
- Mixing tab with extensible shower head (1,50m)
- **3" HMI touchscreen control panel** to control air flow (ON / OFF), lighting, level monitoring for formalin and rinsing system
- **Integrated continuous rinsing system** beneath the work surface, controlled via the touchscreen control panel
- Sink (120 x 120 mm) for **formalin disposal with integrated level sensor and audio-visual alarm displayed via the touchscreen control panel – 10 L canister**
- **Ventilated 25 L plastic waste bin** with left-side swiveling lid

- **Medical rail** for accessories
- **Adjustable feet** to insure table stability
- 4 x **power sockets IP44** with cover
- 1 x **LAN socket**
- Extraction duct connection located opposite the basin, with alternative connection options at the rear of the station (wall connection) or through the underside of the station (floor connection)
- Power supply: **230 V / 50 Hz / 1P / 16 A**
- Air supply: **approx. 800 m<sup>3</sup>/h** | Ø 250 mm | Pressure loss: 50 PA
- Air exhaust: **approx. 1050 m<sup>3</sup>/h** | Ø 280 mm | Pressure loss: 120 PA

## Compatible accessories



Formalin dispensing system with foot pedal

FP000882 | 060.300.001




Water dispensing with foot pedal

FP000661 | 060.300.013




Fixed safety glass side panel for the left or right side of the station

FP001136 | 060.300.023




Document holder

FP000657 | 060.300.002



Keyboard holder and mouse tray (max. 2 kg)

FP001086 | 060.300.003




Monitor holder (max. 11 kg)

FP001085 | 060.300.004



Camera mount with articulated arm and mounting bracket (max. 1,5 kg)

FP001115 | 060.300.019



Splash guard with 3-joint-arm and gooseneck

FP000659 | 060.300.006



Infrared sensor faucet

**FP001093 | 060.300.017**



Waste disposal unit

**FP001092 | 060.300.005**




Magnifying glass with  
LED lightning

**FP001087 | 060.300.008**



Extensible eyes-wash  
shower

**FP001091 | 060.300.009**



Magnetic bar for instruments  
- Length 500 mm

**FP001090 | 060.300.024**



Test tube holder (Ø 20 mm)

**FP001096 | 060.300.016**




Paper dispenser

**FP000658 | 060.300.010**




Gloves dispenser - 1  
size

**FP001118 | 060.300.020**



ORNet Specimen  
Documentation System

**On request**



White PE cutting board  
with non-slip feet

**Material:** Polyethylene

**Dimensions (L/W/H):**  
300 x 400 x 35 mm

**MC002638 | 060.012.011**

**Dimensions (L/I/H):**  
500 x 400 x 35 mm

**MC002637 | 060.012.010**

**Dimensions (L/I/H):**  
600 x 400 x 35 mm

**MC002636 | 060.019.010**




Blue PE cutting board  
with non-slip feet

**Material:** Polyethylene

**Dimensions (L/W/H):**  
500 x 400 x 35 mm

**MC002635 | 060.012.119**



Corian cutting board with groove, marking point and non-slip feet:

**Material:** Corian

**Dimensions (L/W/H):** 500 x 400 x 30 mm

**MC002633 | 060.012.002**



Corian cutting board with groove and non-slip feet

**Material:** Corian

**Dimensions (L/W/H):** 500 x 400 x 30 mm

**MC002634 | 060.012.003**



Cutting board with central black reference point and 30 cm ruler

**Material:** Nabolen 500 (HDPE)

**Dimensions (L/W/H):** 400 x 300 x 25 mm

**FP000655 | 060.300.015**